# visomat

## handy

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Wrist Blood Pressure Monitor

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Instructions for use

**UEBE**  
Germany est.1890

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## Table of contents

<b>A</b>	<b>Intended use</b>	30
<b>B</b>	<b>Safety instructions</b>	
	1. Important patient information	31
	2. Important technical information	32
<b>C</b>	<b>Operating the unit</b>	
	1. Unit description	35
	2. Display	36
	3. Important instructions for use	36
	4. Initial operation of the device	37
	5. Inserting/replacing batteries, battery icon	38
	6. Activating and setting the date/time	38
	7. Attaching the cuff	39
	8. Body posture when taking measurements	40
	9. Measuring blood pressure	41
	10. Classification of the measured values (WHO)	41
	11. Pulse pressure	42
	12. Irregular pulse waves and cardiac arrhythmia	43
	13. Using the memory	43
<b>D</b>	<b>What you should know about blood pressure</b>	
	1. Systolic and diastolic blood pressure values	45
	2. Reasons for measuring different values	45
	3. Why you should measure blood pressure regularly	45
<b>E</b>	<b>Technical information</b>	
	1. Error messages	46
	2. Customer service	47
	3. Technical data	48
	4. Original spare parts and accessories	49
	5. Metrological inspection	49
	6. Explanation of symbols	50
	7. Disposal	50

## Table of contents

<b>F</b>	<b>Maintaining the unit</b>	51
<b>G</b>	<b>Warranty</b>	52

Thank you for choosing the visomat handy wrist blood pressure monitor (also referred to in the following as the unit).

These instructions for use are intended to help the user operate the device safely and efficiently. The unit must be used in accordance with the procedures contained in these instructions for use and must not be used for any other purposes.

It is important that you read all the instructions carefully before using the device. Please pay particular attention to the chapter "Important usage information" on page 36.

This device is intended for persons over 12 years of age for non-invasive self-measurement of the systolic and diastolic blood pressure, determining the pulse rate and calculating the pulse pressure.

This unit uses the oscillometric method for measuring blood pressure and pulse rate.

## 1. Important patient information

- The device may only be used at the wrist and not on any other extremities.
- Blood pressure measurements on children require specialist knowledge! Please consult your doctor if you wish to measure the blood pressure of a child. Do not under any circumstances use the unit on a baby/infant.
- The results of the automatic blood pressure monitor can be distorted by cardiac dysrhythmia, e.g. atrial fibrillation, premature ventricular contraction or arteriosclerosis. Only use the blood measurement device by yourself if you have discussed this with your doctor first.
- During pregnancy and in cases of pre-eclampsia, the device may only be used after consulting your doctor.
- Do not under any circumstances place the cuff over a critical area, e.g. a wound, aneurysm, etc. or on an arm with an arteriovenous shunt. Risk of injury! Any supply via an intravascular access point (infusion) or other medical monitoring devices could possibly be interrupted.
- Do not use the unit without first consulting your doctor if you are undergoing dialysis treatment or taking anticoagulants, platelet aggregation inhibitors or steroids. Internal bleeding may be caused under these conditions.
- The display of the pulse frequency is not suitable for checking the frequency of cardiac pacemakers. Cardiac pacemakers and blood pressure monitors do not influence each other regarding their mode of operation.
- The device is not suitable for use in close proximity to electrosurgical instruments or magnetic resonance imaging equipment.
- Do not use the device in explosive environments such as in the vicinity of flammable anaesthetics or in an oxygen chamber.

- Please wait a few minutes between successive measurements, otherwise the blood flow in the arm is interrupted for too long and injuries could occur.
- The device contains small parts that could be swallowed by children. Do not leave the device unattended with children.
- Store the device out of the reach of children.
- Do not take a blood pressure measurement on the same side as a mastectomy with simultaneous removal of the axilla lymph nodes.
- Discuss the measured values with your doctor. Do not assess the measurement results yourself. **Never change the dose of medicines prescribed by your doctor.**
- Before carrying out your self-measurement, please pay attention to the chapter "Important usage information" on page 36.

## 2. Important technical information

- The unit contains sensitive parts and must be protected from extreme temperature fluctuations, humidity, shocks, dust and direct sunlight.
- Do not knock or drop the device.
- The device is not waterproof. Never immerse it in water or any other liquids. If liquid enters the device, this may result in serious damage and may cause malfunctions.
- Do not apply pressure to the display. Do not place the device with the display facing downwards.
- The unit must only be operated with the cuff and accessories provided (see "Original spare parts and accessories" on page 49). Using different cuffs or accessories can lead to incorrect measurement results.

The warranty will be invalidated if the unit is damaged by third-party accessories!

- A consistently good power supply to the device is necessary for fault-free blood pressure measurements.
  - Only use long-life alkaline batteries (LR03).
  - Always replace all batteries simultaneously during battery replacement. Do not mix new and old batteries or batteries of different types.
  - You will require 2 x 1.5 V batteries. Rechargeable batteries only have a voltage of 1.2 V and are therefore not suitable.
- Never simultaneously touch the device's live electrical parts and the patient.
- If the device is not to be operated for an extended period of time, please remove the batteries. In principle, all batteries can leak.
- Do not remove batteries from the device when the device is switched on.
- Using this device in the vicinity of mobile phones, microwaves or other devices with strong electromagnetic fields can cause malfunctions and inaccurate measurements. During use, maintain a distance of at least 30 cm from such devices.
- Never open or modify the device or the cuff (except for replacing the batteries). If the unit has been opened, it must be subjected to a metrological inspection by an authorised institution.
- To avoid inaccurate readings, please adhere to the intended operating and storage conditions. See "Technical data" on page 48.
- The inflating and measuring procedure can be interrupted by pressing the Start/Stop button or by removing the cuff. The device then stops the inflation procedure and deflates the cuff.


- If the device malfunctions or is faulty, please refer to the troubleshooting guide starting on page 46, or contact customer service (see page 47).

## 1. Unit description

- 1 Display
- 2 Start/Stop button
- 3 Memory button
- 4 Battery compartment (Removable battery compartment cover)

## 2. Display

- |                                                                        |                                       |
|------------------------------------------------------------------------|---------------------------------------|
| <b>1</b> SYS = Systole                                                 | <b>5</b> Battery check display        |
| <b>2</b> DIA = Diastole                                                | <b>6</b> Inflating / Deflating        |
| <b>3</b> PUL 1/min = Pulse<br>Calculated pulse frequency per<br>minute | <b>7</b> WHO classification (Page 41) |
| <b>4</b> Pulse signal display or irregu-<br>lar pulse waves            | <b>8</b> Memory identifier            |
| <b>9</b> Storage location or date/<br>time                             |                                       |

## 3. Important instructions for use

The measurement results of automated blood pressure monitors can be influenced by the measuring location, posture, previous exertion and general physical condition. Please pay attention to the usage information in order to obtain accurate readings.

- Do not consume alcohol, caffeine or smoke for at least one hour before measuring.

- Rest for at least 5 minutes before taking the measurement. Depending on the previous degree of exertion, you may even need to rest for up to an hour.
- Expose the wrist. Under no circumstances may clothes hinder the blood flow to or from the arm, as this affects the blood pressure at the measuring point and can lead to incorrect readings.
- Do not move or talk during measuring.
- Breathe slowly and deeply. Do not hold your breath.
- Check the display for irregular pulse signals. If required, repeat the measurement when conditions have improved.
- In order to exclude side differences and to obtain comparable measurement results, it is important to always take the measurement on the same arm. Ask your doctor which side it is advisable for you to measure.
- Always take your daily blood pressure measurements at around the same time.
- Blood pressure is not a fixed value. It can go up or down by more than 20 mmHg in patients within a few minutes.

## 4. Initial operation of the device

Insert the batteries provided into the device.

The device has a date/time function that can be activated as required (see page 38). This function is deactivated on delivery.

When the date/time function is activated, the measurement results are saved along with their date and time. This can make it easier to view the readings at a later point.

## 5. Inserting/replacing batteries, battery icon

- Remove the cover of the battery compartment on the side of the unit.
- Remove the old batteries from the device and insert the new ones. Please ensure correct polarity (marking in battery compartment).
- Close the battery compartment by clipping the battery cover back into the device.

When the “empty battery ” symbol first appears in the display, you are still able to take around 30 measurements. Please replace the batteries within this time frame.

## 6. Activating and setting the date/time

### Activating and deactivating the date/time function

To activate the date/time function, press the Start/Stop button and then immediately, while the full display is visible, the Memory button. This function can be deactivated at any time in the same way.

When the date/time function is activated, the device continues to display the time even when it is switched off. The power consumption for this is extremely low.

### Setting the date/time

The four-digit number representing the year appears in the display first (Figure 1). Press the Memory button to set the year. Press the Start/Stop button to confirm the selected year and proceed to setting the date (Figure 2).

Figure 1

Set the month with the Memory button and confirm this with the Start/Stop button. Set the day, hours and minutes in the same way (Figure 2, Figure 3).

Figure 2

Figure 3

### Changing the date/time


To change the date and time, take the batteries out and wait until the display goes off. Then re-insert the batteries and reset the date and time.

After replacing the batteries, you will have to reset the date and time.

## 7. Attaching the cuff

The special double-preformed comfort cuff is easy to attach and ensures a perfect fit during measurement. The materials used are tested for skin compatibility according to ISO 10993.

- Clear the wrist from all clothing or watches/jewellery since this may disturb the blood flow in the arm as it may impair the blood pressure at the place of measurement.
- Slide the device onto the centre of the inside of the forearm, directly on the wrist. The palm of your hand should be facing upwards. The distance between the wrist and the upper edge of the cuff should be 5–10 mm.
- Fasten the cuff using the hook and loop fastening.

The shape of the cuff is optimized to sit very tightly in order to avoid possible measurement errors.

WHO 2003	Systolic pressure = Upper value mmHg**	Diastolic pressure = Lower value mmHg**
Red	from 140*	from 90*
Yellow	120 to 139	80 to 89
Green	below 120	below 80

\* It is sufficient for one of the values to be raised.

\*\* millimetres of mercury

The categorisation according to WHO is stored with the readings and can be called up again together with the readings from the memory.

Depending on age, weight and general health, the blood pressure values can be different. Only a doctor can determine the right blood pressure range for you and assess whether your blood pressure has reached a dangerous level for you. Discuss your blood pressure values with your doctor. **Never change the dose of medicines prescribed by your doctor!**

## 11. Pulse pressure

The pulse pressure – not to be confused with the pulse beat – provides an indication of the elasticity of the blood vessels. A stiff system of vessels can have a negative effect on the cardiovascular system. Studies show that the cardiovascular risk increases if the pulse pressure value is permanently higher than 65 mmHg.

The heart works in two phases, the contraction phase (systole) and the relaxation phase (diastole). The pressure difference between systole and diastole is called the pulse pressure or pulse amplitude. The higher the pulse pressure, the harder the vessels.

High pulse pressure	above 65 mmHg
Increased pulse pressure	55 to 65 mmHg
Normal pulse pressure	below 55 mmHg

If your pulse pressure is permanently above 65 mmHg, consult your doctor.

## 12. Irregular pulse waves and cardiac arrhythmia

If the device detects irregular pulses during the measurement, the ♥ symbol is displayed once the measurement is complete. This can be set off by an irregular heartbeat (arrhythmia), disturbances caused by movement, talking or even breathing in deeply. The symbol is stored with the respective measurement.

If the symbol appears more frequently, this may be caused by an irregular heartbeat and must be discussed with your doctor!

Measurement results accompanied by irregular pulses should therefore be viewed as critical and repeated under more favourable conditions.

**Caution:** Do not confuse the irregular pulse display ♥ with the flashing pulse display ♥ that indicates the heart rate during every measurement (see page 41).

## 13. Using the memory

Measured results are automatically stored in the memory. The memory can store up to 60 results and the average value.

When more than 60 measured values have been stored, the oldest value (No. 60) is deleted to allow the latest value (No. 1) to be recorded.

### Data retrieval

Press the Memory button to call up data. The average value of the stored results is displayed with "A", the values for systole, diastole and pulse and the pulse pressure (PP) are shown alternately in the display (Figure 1, Figure 2).


Figure 1

Figure 2

When the Memory button is pressed again, the value measured last (No. 1) appears. When the date/time function is activated, systole, diastole and pulse are displayed alternately with memory number (Figure 3) and date (Figure 4), followed by the pulse pressure (PP) and the time (Figure 5).

Figure 3

Figure 4

Figure 5

Press the Memory button repeatedly to call up further measured values.

The stored data remain visible in the display for around 30 seconds. The unit then switches itself off.

### Deleting data

To delete individual measured values, press the Memory button repeatedly to call up the desired measured value.

Then press and hold the Memory button again for 8–10 seconds until the measured value flashes several times and then goes out.

To delete the entire memory, retrieve the average value then press and hold the Memory button again for 8–10 seconds until the measured value flashes several times and then goes out.

If no value is stored, there is no display.

## 1. Systolic and diastolic blood pressure values

The cardiovascular system has the important function of supplying all organs and tissues in the body with sufficient amounts of blood and of transporting metabolites. For this, the heart contracts and expands at a regular rate of about 60 to 80 times per minute. The pressure of the flowing blood on the artery walls caused by the heart contracting is termed systolic. The pressure in the ensuing relaxation phase, when the heart re-fills with blood, is termed diastolic. During daily measurement you determine both values.

## 2. Reasons for measuring different values

Our blood pressure responds to internal and external influences like a sensitive measuring instrument. Affected by mental, physical and environmental influences, it varies all the time and never remains constant. Reasons for fluctuating blood pressure values can be: moving, speaking, eating, consuming alcohol or nicotine, nervousness, inner tension, emotions, room temperature, recent urination or bowel movements, environmental influences such as movements and sounds, etc. Even changes in the weather and climate can have an effect on your blood pressure. This explains why values measured at the doctor are often higher than those you obtain at home in your usual environment.

## 3. Why you should measure blood pressure regularly

Even the time of day has an effect on your blood pressure. The values are at their highest during the day. In the course of the afternoon and in the evening, they drop slightly. They are low while you are sleeping, but rise again relatively quickly once you get up. One-off and irregular measurements therefore say little about your actual blood pressure. A reliable assessment is possible only when measurements are taken regularly. Discuss the measurement values with your doctor.

## 1. Error messages

Failure encountered	Possible cause	Corrective action
Display Err - 300	Excessive cuff pressure. Cuff inflated to maximum due to movement of the arm or body.	-Repeat the measurement -Do not move your arm -Do not speak
Display Err - 1	Measurement error; measurement could not be carried out due to movement or speaking.	-Repeat the measurement -Do not move your arm -Do not speak
Display Err - 2	Inflation error, pressure built up too quickly or too slowly. Cuff is not correctly fitted.	-Repeat the measurement -Do not move your arm -Do not speak -Check wrist circumference.
Display Err - 3	Deflation rate too high or too low. Movement or talking during measurement.	-Repeat the measurement -Do not move your arm -Do not speak
Display 0 0	Start/Stop button accidentally pressed while batteries were being changed.	Switch unit off and then on again with Start/Stop button.
	Irregular heartbeat, disruptive movements, shaking, wobbling (objects), breathing in deeply, etc.	Repeat measurement after 3-5 minutes rest. Possible influences on the measurement results (depending on the severity of the arrhythmia) must be discussed with your doctor!
Measured values/ pulse values are extremely high or low.	Too little rest before the measurement. Movement or talking during measurement.	Repeat the measurement after a break of 3-5 minutes. Pay attention to the usage information on page 36.
	Incorrect measurement posture	Repeat the measurement and maintain the correct posture. (Page 40)
Consecutive measurements produce different measured values (readings).	Blood pressure is not a fixed quantity. It may vary in the upward or downward direction within a few minutes by more than 20 mmHg.	In order to obtain comparable results, always carry out the measurements under the same conditions.

Failure encountered	Possible cause	Corrective action
Measurement was interrupted, the cuff deflated to 0 and the measurement restarted.	Device detects a fault or that the diastolic value is too low. Movement during measurement	Device shows no fault. Repeat measurement If necessary, interrupt measurement, measure again after a break of 5 minutes.
No display after unit is switched on	Batteries inserted incorrectly.	Check position of batteries.
	Batteries flat. Battery contacts dirty.	Change batteries. Clean battery contacts with a dry cloth.
Nothing appears in the display after pressing the Memory button	No measured values have been stored.	Repeat measurement
Measurement interrupted	Batteries flat.	Change batteries.
The  symbol appears during measurement but the measurement continues.	Low batteries	You can still perform approx. 30 measurements, but should keep spare batteries available.
Date/time not shown in display	Function accidentally deactivated	Activate date and time as described on page 38.

## 2. Customer service

Device repairs may only be performed by the manufacturer or an expressly authorised body. Please contact:

UEBE Medical GmbH  
 Bürgermeister-Kuhn-Straße 22  
 97900 Kulsheim, Germany  
 info@uebe.com  
 www.uebe.com

### 3. Technical data

Model:	visomat handy
Size:	approx. W = 70 mm x H = 70 mm x D = 80 mm with cuff
Weight:	121 g with cuff, not including batteries
Display:	LCD display (liquid crystal display) for measured values and check displays
Memory:	60 measured values (stored automatically) and average value (A)
Measurement procedure:	Oscillometric determination of systole, diastole and pulse
Reference procedure of clinical testing:	Auscultatory measurement
Inflating pressure:	Automatically controlled
Pressure display range:	0-300 mmHg
Measurement range:	Systolic: 50-250 mmHg Diastolic: 40-180 mmHg Pulse measurement: 40-160 pulses/min
Measurement precision:	Blood pressure measurement (cuff pressure): $\pm 3$ mmHg Pulse rate: $\pm 5\%$
Power supply:	Battery type: 2 x 1.5 V, AAA, LR 03, durability: over 700 measurements in 2 years
Power consumption:	2 W
Cuff:	Cuff type visomat H5, wrist circumference 12.5 cm - 21.5 cm
Operating conditions:	Ambient temperature 10°C to 40°C, relative humidity 15% to 85% non-condensing, air pressure 700 to 1060 hPa
Storage and transportation conditions:	Ambient temperature -20 °C to 60 °C, relative humidity 15 to 95%, non-condensing
Time to adjust from minimum/maximum storage temperature to operating temperature:	2 hours
Inflation and deflation:	Regulated automatically
Protection class IP:	IP 20: Protected against solid foreign particles with a diameter of more than 12.5 mm, no protection against water.

Protection against electric shock:	Internal power supply, applied part type BF (cuff)
Expected service life:	5 years
Operating mode:	Continuous operation
Classification:	Internally powered ME Equipment

### 4. Original spare parts and accessories

The following original spare parts and accessories are available from specialist dealers:

- Cuff type H5 (for wrist circumference 12.5-21.5 cm)  
Part no. 2105001, PZN-06414493

Subject to technical modifications.

### 5. Metrological inspection

Generally, it is recommended that a metrological inspection is performed every two years. However, professional users in Germany are obligated to do so in accordance with the "Medical Products Operator Ordinance" (Medizinprodukte-Betreiberverordnung). This can be performed either by UEBE Medical GmbH, an authority responsible for metrology or an authorised maintenance service. Please refer to your national regulations.

Upon request, responsible authorities or authorised maintenance services receive a "Test instruction for metrological inspection" from the manufacturer.

Please only submit the device for metrological inspection together with the cuff and the instructions for use.

Important: No modifications, e.g. opening the device (except to replace the batteries), may be made to this device without the manufacturer's permission.

## 6. Explanation of symbols

This product complies with Council Directive 93/42/EEC of 5 September 2007 concerning medical devices and bears the mark CE 0123 (TÜV SÜD Product Service GmbH).

Degree of protection against electric shock TYPE BF

Consult the instructions for use

Keep dry

Device serial number

Reference number = item number

**IP20**

Protected against solid foreign particles with a diameter of more than 12.5 mm, no protection against water.

Manufacturer

## 7. Disposal

Batteries and technical appliances must not be disposed of with domestic waste, but should be handed in at the appropriate collection and disposal points.

- Keep the device clean. Check for cleanliness after use. Please use a soft, dry cloth for cleaning. Do not use benzene, thinners or other strong solvents.
- The cuff may absorb sweat and other fluids. Inspect the cuff for marks or discolourations after each use. Please use a synthetic detergent for cleaning and carefully wipe the surface. Do not scrub or machine-wash. Air-dry thoroughly.
- Make sure when storing the unit that no heavy objects are left lying on the unit.
- Do not apply force when folding or bending the cuff.
- If the device has been stored at temperatures below 10°C or above 40°C, allow it to stand at room temperature for at least 2 hours before using it.
- If the device is not to be operated for an extended period of time, please remove the batteries. Store the batteries out of the reach of children.

The device has been manufactured and tested with all due care. Nevertheless, in the event of defects upon delivery, we provide a warranty under the following conditions:

During the warranty period of 5 years from the date of purchase, we remedy such defects at our discretion and at our expense in our plant through repair work or replacement delivery of a defect-free device.

The warranty does not cover normal wear and tear of wearing parts or damage caused by failure to observe the instructions for use, improper handling (e.g. unsuitable power sources, breakage, leaking batteries) and/or disassembly of the device by the buyer. Furthermore, the warranty does not constitute grounds for asserting claims for damages against us.

Warranty claims can only be asserted during the warranty period and upon presentation of the proof of purchase. In the case of a warranty claim, the device along with the proof of purchase and a description of the complaint must be sent to:

UEBE Medical GmbH  
Service-Center  
Zum Schlag 18  
97877 Wertheim, Germany

The cost of returning the device to our factory shall be borne by the sender. Complaints that are sent back without prepayment will not be accepted by UEBE.

The statutory claims and rights of the buyer against the seller (for example, defect claims, producer liability) are not restricted by this warranty.

**Please note: In the case of a warranty claim, please make sure to enclose the proof of purchase.**

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© Copyright 2020 UEBE Medical GmbH